HIPRA

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Practical approach to
intradermal vaccination, Hipradermic

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Alex WU, Regional Technical and Marketing Manager Asia & Oceania

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It's time to focus on what really matters

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快適さ & 正確さ Comfort&Precise

皮内無針投与器 The intradermal needle-free device.


We need to replace this video.
Could you put on Piglet's video?

トリガーなし: 使用者の疲労感を軽減し、投与時間も短縮
Automatic injection without a trigger,
so reduces user muscle fatigue and administration time.

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Hands-free option available for an even easier piglet's vaccination.

Providing prevention
with innovative solutions
for livestock

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
On the volume injected: **0.2 ml**

On the dispersion profile and tissue depth

一定の投与量を保つ Consistency with the does administered

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
利便性 Convenience

ユーザー想いのデバイス：
User-friendly device for carrying out
vaccinations easily and quickly.

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The intradermal injection with no needle is a **less invasive technique**,
which means **better performance**

子豚は痛いとき高周波の声で鳴く
HIGH FREQUENCIES ARE LINKED TO PAIN EPISODES (PUPPE ET AL., 2005)

高周波High frequencies (> 1000 HZ)

低周波Low frequencies (< 1000 HZ)


Intramuscular vs. Intradermic Needle-Free Vaccination in Piglets: Relevance for Animal Welfare Based on an Aversion Learning Test and Vocalizations

OPEN ACCESS

Edited by:
Kathleen Mary O'Driscoll,

Antoni Dalmau^{1*}, Almudena Sánchez-Matamoros², Jorge M. Molina², Aida Xercavins¹, Aranzazu Varvaró-Porter¹, Israel Muñoz¹, Xènia Moles¹, Berta Baulida¹, Emma Fàbrega¹, Antonio Velarde¹, Joaquim Palliserà¹, Alba Puigredon² and Alexandra Contreras-Jodar¹

¹ Animal Welfare Program, Institute of Agrifood Research and Technology, Barcelona, Spain; ² HIPRA, Girona, Spain

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The intradermal injection with no needle is a **less invasive technique**,
which means **better performance**

Intradermal needle-free vaccination in piglets: effects on welfare and nursery performance

Daniela Bruna Ferrandin¹; Tatiana C. G. Dutra de Souza¹; Gabriela Ibanez Barbosa¹; Gabriel Peixotto Faria¹; Mickael Javan Camello¹; Ramon Jorda Casadevall²; Lorena Nodar Balseiro²

¹HIPRA Saúde Animal Brasil, Porto Alegre, RS ²Laboratorios HIPRA, Girona, Spain., daniela.ferrandin@hipra.com

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*Significant differences between Intradermal group and Intramuscular group were observed (Linear Regression test, $p < 0.001$; was performed)

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
Better CMI response and lower immunosuppression

With a **minimal antigen dose (0,2 ml)** we get an
improved immune response

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Adthakorn Madapong et., al. 2020 *Vet Microbiol.* 244:108655

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
No needles are involved so **the risk of disease**

transmission between animals is reduced and internal biosecurity is improved

Evaluation of ASF and PRRS virus transmission between pigs when using conventional needles and a needle-free device

Dachrit Nilubol^{1,2,*}, Joel Miranda³, Salvador Romero³, Sittikorn Traiyarach³, Angkana Tantituvanont^{2,4}, Dante Palabrica⁵

¹Swine Viral Evolution and Vaccine Development Research Unit and ²Department of Veterinary Microbiology, Faculty of Veterinary Science, Chulalongkorn University, Thailand, ³HIPRA, Amer (Girona), Spain, ⁴Department of Pharmaceutical and Industrial Pharmacies, Faculty of Pharmaceutical Sciences, Chulalongkorn University, Thailand, ⁵Robina Farms Diagnostic Laboratory, Universal Corn Products Compound, Philippines,

*Corresponding author: dachrit@gmail.com

	GROUP	SEROCONVERSION OF SENTINELS TO ASF AND PRRS AFTER ID AND IM INOCULATION 感染実験後に抗体陽転を示した感受性豚の頭数				
		0 DPI	7 DPI	14 DPI	21 DPI	28 DPI
ASF	Needle 注射針	0/6 (Negative)	0/6 (Negative)	4/6 (Positive)	6/6 (Positive)	6/6 (Positive)
	Hipradermic®	0/6 (Negative)	0/6 (Negative)	0/6 (Negative)	0/6 (Negative)	0/6 (Negative)
PRRS	Needle 注射針	0/6 (Negative)	0/6 (Negative)	2/6 (Positive)	6/6 (Positive)	6/6 (Positive)
	Hipradermic®	0/6 (Negative)	0/6 (Negative)	0/6 (Negative)	0/6 (Negative)	0/6 (Negative)

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